

Product datasheet for **GA210430**

Mouse Ttpal activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ttpal CRISPRa kit - CRISPR gene activation of mouse tocopherol (alpha) transfer protein-like
Symbol:	Ttpal
Components:	GA210430G1, Ttpal gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210430G2, Ttpal gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210430G3, Ttpal gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D3D0</u>
Synonyms:	3110080A02Rik; 5830472M02Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)