

Product datasheet for **GA210446**

Mouse Lamp5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lamp5 CRISPRa kit - CRISPR gene activation of mouse lysosomal-associated membrane protein family, member 5
Symbol:	Lamp5
UniProt ID:	Q9D387
Synonyms:	3110035N03Rik; 6330527O06Rik; BAD-LAMP; LAMP-5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)