

Product datasheet for **GA210470**

Mouse Tsen54 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tsen54 CRISPRa kit - CRISPR gene activation of mouse tRNA splicing endonuclease subunit 54
Symbol:	Tsen54
UniProt ID:	Q8C2A2
Synonyms:	0610034P02Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)