


Product datasheet for **GA210501**

Mouse Ift43 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ift43 CRISPRa kit - CRISPR gene activation of mouse intraflagellar transport 43
Symbol:	Ift43
Components:	GA210501G1, Ift43 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210501G2, Ift43 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210501G3, Ift43 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9DA69
Synonyms:	1700019E19Rik; R75064
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)