

Product datasheet for GA210519

Mouse Haus8 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Haus8 CRISPRa kit - CRISPR gene activation of mouse 4HAUS augmin-like complex, subunit 8
Symbol:	Haus8
Components:	GA210519G1, Haus8 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210519G2, Haus8 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210519G3, Haus8 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q99L00</u>
Synonyms:	2410004L22Rik; Hice1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

