


Product datasheet for **GA210569**

Mouse Wdr38 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr38 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 38
Symbol:	Wdr38
Components:	GA210569G1, Wdr38 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210569G2, Wdr38 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210569G3, Wdr38 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D994</u>
Synonyms:	1700123D08Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)