

Product datasheet for GA210594

Mouse Efr3a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Efr3a CRISPRa kit - CRISPR gene activation of mouse EFR3 homolog A
Symbol:	Efr3a
Components:	GA210594G1, Efr3a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210594G2, Efr3a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210594G3, Efr3a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BG67</u>
Synonyms:	A130089M23Rik; BB071175; C76891; C920006C10Rik; D030063F01Rik; mKIAA0143
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

