

Product datasheet for **GA210615**

Mouse 2410137M14Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	2410137M14Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 2410137M14 gene
Symbol:	2410137M14Rik
Components:	GA210615G1, 2410137M14Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210615G2, 2410137M14Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210615G3, 2410137M14Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	RP23-98F21.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)