

Product datasheet for **GA210617**

Mouse Tmem234 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem234 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 234
Symbol:	Tmem234
Components:	GA210617G1, Tmem234 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210617G2, Tmem234 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210617G3, Tmem234 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8R1E7</u>
Synonyms:	1500002D11Rik; 2510006D16Rik; 4933407D05Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)