

## Product datasheet for **GA210656**

### Mouse Spata46 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | Spata46 CRISPRa kit - CRISPR gene activation of mouse spermatogenesis associated 46                                                                                                                                                                                                                              |
| Symbol:       | Spata46                                                                                                                                                                                                                                                                                                          |
| Components:   | GA210656G1, Spata46 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA210656G2, Spata46 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA210656G3, Spata46 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q4FZF2</a></u>                                                                                                                                                                                                                                                                                    |
| Synonyms:     | 1700015E13Rik                                                                                                                                                                                                                                                                                                    |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)