

Product datasheet for **GA210657**

Mouse Tsacc activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tsacc CRISPRa kit - CRISPR gene activation of mouse TSSK6 activating co-chaperone
Symbol:	Tsacc
Components:	GA210657G1, Tsacc gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210657G2, Tsacc gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210657G3, Tsacc gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9DA44</u>
Synonyms:	1700021C14Rik; SIP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)