


## Product datasheet for **GA210659**

### Mouse Ifi2712a activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                     |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                        |
| Description:  | Ifi2712a CRISPRa kit - CRISPR gene activation of mouse interferon, alpha-inducible protein 27 like 2A                                                                                                                                                                                                               |
| Symbol:       | Ifi2712a                                                                                                                                                                                                                                                                                                            |
| Components:   | GA210659G1, Ifi2712a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA210659G2, Ifi2712a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA210659G3, Ifi2712a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q8R412</a>                                                                                                                                                                                                                                                                                              |
| Synonyms:     | 2310061N23Rik; Ifi27; lsg12; lsg12(b1)                                                                                                                                                                                                                                                                              |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                             |

#### Product images:



[View online »](#)