

Product datasheet for **GA210683**

Mouse Ticrr activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ticrr CRISPRa kit - CRISPR gene activation of mouse TOPBP1-interacting checkpoint and replication regulator
Symbol:	Ticrr
UniProt ID:	Q8BQ33
Synonyms:	5730590G19Rik; BC028450; Lrrk1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)