

Product datasheet for **GA210688**

Mouse Tstd3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tstd3 CRISPRa kit - CRISPR gene activation of mouse thiosulfate sulfurtransferase (rhodanese)-like domain containing 3
Symbol:	Tstd3
Components:	GA210688G1, Tstd3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210688G2, Tstd3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210688G3, Tstd3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9D0B5
Synonyms:	2610029I01Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)