

Product datasheet for **GA210707**

Mouse Tanc2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tanc2 CRISPRa kit - CRISPR gene activation of mouse tetratricopeptide repeat, ankyrin repeat and coiled-coil containing 2
Symbol:	Tanc2
Components:	GA210707G1, Tanc2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210707G2, Tanc2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210707G3, Tanc2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A2A690
Synonyms:	3526402J09Rik; 5730590C14Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)