

Product datasheet for **GA210709**

Mouse Gbp111 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Gbp111 CRISPRa kit - CRISPR gene activation of mouse GC-rich promoter binding protein 1-like 1
Symbol:	Gbp111
Components:	GA210709G1, Gbp111 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210709G2, Gbp111 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210709G3, Gbp111 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6NZP2
Synonyms:	5330440M15Rik; BC002292
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)