


## Product datasheet for **GA210771**

### Mouse Tcf1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | Tcf1 CRISPRa kit - CRISPR gene activation of mouse TRPM8 channel-associated factor 1                                                                                                                                                                                                                    |
| Symbol:       | Tcf1                                                                                                                                                                                                                                                                                                    |
| Components:   | GA210771G1, Tcf1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA210771G2, Tcf1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA210771G3, Tcf1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u>Q8BNE1</u>                                                                                                                                                                                                                                                                                           |
| Synonyms:     | 2810407D09Rik; 3321401G04Rik; A230020K05Rik; AI848529; Fam115a                                                                                                                                                                                                                                          |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)