

Product datasheet for GA210784

Mouse Odad3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ccdc151 CRISPRa kit - CRISPR gene activation of mouse coiled-coil domain containing 151
Symbol:	Odad3
Components:	GA210784G1, Odad3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210784G2, Odad3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210784G3, Odad3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BSN3</u>
Synonyms:	AI64444I5; b2b1885Clo; C330001K17Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

