

Product datasheet for **GA210793**

Mouse Snapc3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Snapc3 CRISPRa kit - CRISPR gene activation of mouse small nuclear RNA activating complex, polypeptide 3
Symbol:	Snapc3
Components:	GA210793G1, Snapc3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210793G2, Snapc3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210793G3, Snapc3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9D2C9
Synonyms:	1810020H02Rik; 4930558A07Rik; 5031401C21Rik; A1414457; E030018J20Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)