


## Product datasheet for **GA210794**

### Mouse C330007P06Rik activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | C330007P06Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA C330007P06 gene   |
| Symbol:       | C330007P06Rik  |
| Components:   | GA210794G1, C330007P06Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA210794G2, C330007P06Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA210794G3, C330007P06Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q8VDP2</a></u>  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)  |

#### Product images:



[View online »](#)