

Product datasheet for **GA210798**

Mouse Arhgef38 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Arhgef38 CRISPRa kit - CRISPR gene activation of mouse Rho guanine nucleotide exchange factor (GEF) 38
Symbol:	Arhgef38
Components:	GA210798G1, Arhgef38 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210798G2, Arhgef38 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210798G3, Arhgef38 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80VK6</u>
Synonyms:	9130221D24Rik; D630013G24Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)