

Product datasheet for **GA210816**

Mouse Elp4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Elp4 CRISPRa kit - CRISPR gene activation of mouse elongator acetyltransferase complex subunit 4
Symbol:	Elp4
Components:	GA210816G1, Elp4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210816G2, Elp4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210816G3, Elp4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9ER73</u>
Synonyms:	A330107A17Rik; Paxneb
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)