

Product datasheet for GA210817

Mouse Ermn activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ermn CRISPRa kit - CRISPR gene activation of mouse ermin, ERM-like protein
Symbol:	Ermn
Components:	GA210817G1, Ermn gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210817G2, Ermn gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210817G3, Ermn gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5EBJ4
Synonyms:	A330104H05Rik; A1854460; ermin; mKIAA1189
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

