

Product datasheet for **GA210818**

Mouse Csrnp3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Csrnp3 CRISPRa kit - CRISPR gene activation of mouse cysteine-serine-rich nuclear protein 3
Symbol:	Csrnp3
Components:	GA210818G1, Csrnp3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210818G2, Csrnp3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210818G3, Csrnp3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P59055</u>
Synonyms:	A330102K23Rik; Csrnp3; CSRNP-3; Mbu1; taip-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)