

Product datasheet for **GA210827**

Mouse *Lrrc42* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Lrrc42</i> CRISPRa kit - CRISPR gene activation of mouse leucine rich repeat containing 42
Symbol:	<i>Lrrc42</i>
Components:	GA210827G1, <i>Lrrc42</i> gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210827G2, <i>Lrrc42</i> gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210827G3, <i>Lrrc42</i> gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8R2U7
Synonyms:	A930011F22Rik; AA536996
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)