

Product datasheet for **GA210850**

Mouse Rpgrip1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Rpgrip1 CRISPRa kit - CRISPR gene activation of mouse retinitis pigmentosa GTPase regulator interacting protein 1
Symbol:	Rpgrip1
Components:	GA210850G1, Rpgrip1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210850G2, Rpgrip1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210850G3, Rpgrip1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9EPQ2</u>
Synonyms:	0610005A07Rik; 4930401L23Rik; 4930505G06Rik; A930002K18Rik; AA415034; nmf247
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)