

Product datasheet for **GA210892**

Mouse Lpar4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lpar4 CRISPRa kit - CRISPR gene activation of mouse lysophosphatidic acid receptor 4
Symbol:	Lpar4
UniProt ID:	<u>Q8BLG2</u>
Synonyms:	5730485F04Rik; Gpr23; LPA4; p2y9
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)