

Product datasheet for **GA210934**

Mouse H2bu2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hist3h2ba CRISPRa kit - CRISPR gene activation of mouse histone cluster 3, H2ba
Symbol:	H2bu2
Components:	GA210934G1, H2bu2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210934G2, H2bu2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210934G3, H2bu2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D2U9</u>
Synonyms:	1500011O09Rik; AI413321
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)