


Product datasheet for **GA210935**

Mouse Gpr108 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Gpr108 CRISPRa kit - CRISPR gene activation of mouse G protein-coupled receptor 108
Symbol:	Gpr108
Components:	GA210935G1, Gpr108 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210935G2, Gpr108 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210935G3, Gpr108 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q91WD0</u>
Synonyms:	1810015L19Rik; AA589464; C79132; Lustr2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)