

Product datasheet for **GA210948**

Mouse *Icam4* activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | <i>Icam4</i> CRISPRa kit - CRISPR gene activation of mouse intercellular adhesion molecule 4, Landsteiner-Wiener blood group |
| Symbol: | <i>Icam4</i> |
| UniProt ID: | Q9ERM2 |
| Synonyms: | 1810015M19Rik; Cd242; ICAM-4 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)