

Product datasheet for GA210949

Mouse Snrnp25 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Snrnp25 CRISPRa kit - CRISPR gene activation of mouse small nuclear ribonucleoprotein 25 (U11/U12)
Symbol:	Snrnp25
Components:	GA210949G1, Snrnp25 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210949G2, Snrnp25 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210949G3, Snrnp25 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8VIK1</u>
Synonyms:	3300001G02Rik; AL033324; C16ORF33
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



View online »