

## Product datasheet for **GA210984**

### Mouse Utp23 activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | Utp23 CRISPRa kit - CRISPR gene activation of mouse UTP23 small subunit processome component |
| Symbol:       | Utp23  |
| UniProt ID:   | <a href="#">Q9CX11</a>   |
| Synonyms:     | 1700010I21Rik; AI662478; BB238373; D530033C11Rik   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                      |

#### Product images:



[View online >](#)