

Product datasheet for **GA210989**

Mouse Uvr^{ag} activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Uvr ^{ag} CRISPRa kit - CRISPR gene activation of mouse UV radiation resistance associated gene
Symbol:	Uvr ^{ag}
UniProt ID:	Q8K245
Synonyms:	9530039D02Rik; AI648770; BB124205; Uvr ^{ag} 1; Uvr ^{ag} l
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)