

Product datasheet for **GA211030**

Mouse Tepsin activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tepsin CRISPRa kit - CRISPR gene activation of mouse TEPSIN, adaptor related protein complex 4 accessory protein
Symbol:	Tepsin
Components:	GA211030G1, Tepsin gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211030G2, Tepsin gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211030G3, Tepsin gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3U3N6</u>
Synonyms:	2410002I01Rik; Enthd2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)