

Product datasheet for GA211044

Mouse Ttc30a1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ttc30a1 CRISPRa kit - CRISPR gene activation of mouse tetratricopeptide repeat domain 30A1
Symbol:	Ttc30a1
Components:	GA211044G1, Ttc30a1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211044G2, Ttc30a1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211044G3, Ttc30a1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q99J38</u>
Synonyms:	4930506L13Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

