

Product datasheet for **GA211049**

Mouse P2ry10 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2ry10 CRISPRa kit - CRISPR gene activation of mouse purinergic receptor P2Y, G-protein coupled 10
Symbol:	P2ry10
UniProt ID:	Q8BFU7
Synonyms:	5830408N17Rik; Lypsr2; P2Y10
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)