


Product datasheet for **GA211102**

Mouse Sh3bp5l activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sh3bp5l CRISPRa kit - CRISPR gene activation of mouse SH3 binding domain protein 5 like
Symbol:	Sh3bp5l
Components:	GA211102G1, Sh3bp5l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211102G2, Sh3bp5l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211102G3, Sh3bp5l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q99LH9
Synonyms:	2310074E09Rik; BC003251; mKIAA1720
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)