

Product datasheet for **GA211114**

Mouse Sptbn4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sptbn4 CRISPRa kit - CRISPR gene activation of mouse spectrin beta, non-erythrocytic 4
Symbol:	Sptbn4
Components:	GA211114G1, Sptbn4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211114G2, Sptbn4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211114G3, Sptbn4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	1700022P15Rik; 5830426A08Rik; dyn; lnd; nmf261; nmf379; qv; ROSA62; SpbIV; Spnb4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)