

Product datasheet for GA211115

Mouse Kcnip4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Kcnip4 CRISPRa kit - CRISPR gene activation of mouse Kv channel interacting protein 4
Symbol:	Kcnip4
UniProt ID:	Q6PHZ8
Synonyms:	AV032399; Calp; Calp250; KChIP4; KChIP4a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)