

## Product datasheet for **GA211139**

### Mouse H1f1 activation kit by CRISPRa

#### Product data:

<b>Product Type:</b>	CRISPRa Kits
<b>Description:</b>	Hist1h1a CRISPRa kit - CRISPR gene activation of mouse histone cluster 1, H1a
<b>Symbol:</b>	H1f1
<b>Components:</b>	GA211139G1, H1f1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211139G2, H1f1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211139G3, H1f1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
<b>UniProt ID:</b>	<u><a href="#">P43275</a></u>
<b>Synonyms:</b>	H1.1; H1a; H1f1; H1var3
<b>Format:</b>	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

