

Product datasheet for **GA211144**

Mouse Slc16a3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc16a3 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 16 (monocarboxylic acid transporters), member 3
Symbol:	Slc16a3
Components:	GA211144G1, Slc16a3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211144G2, Slc16a3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211144G3, Slc16a3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	P57787
Synonyms:	Mct3; Mct4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)