

Product datasheet for **GA211145**

Mouse Kank3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Kank3 CRISPRa kit - CRISPR gene activation of mouse KN motif and ankyrin repeat domains 3
Symbol:	Kank3
Components:	GA211145G1, Kank3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211145G2, Kank3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211145G3, Kank3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Z1P7</u>
Synonyms:	0610013D04Rik; Ankrd47; D17Ert288e; NG28
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)