

Product datasheet for GA211150

Mouse Tlnrd1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tlnrd1 CRISPRa kit – CRISPR gene activation of mouse talin rod domain containing 1
Symbol:	Tlnrd1
Components:	GA211150G1, Tlnrd1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211150G2, Tlnrd1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211150G3, Tlnrd1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9ERE8
Synonyms:	AW061151
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

