

Product datasheet for **GA211158**

Mouse Dtx3 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Dtx3 CRISPRa kit - CRISPR gene activation of mouse deltex 3, E3 ubiquitin ligase |
| Symbol: | Dtx3 |
| UniProt ID: | Q80V91 |
| Synonyms: | Deltex3; mg310 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)