

Product datasheet for **GA211169**

Mouse Arl4d activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Arl4d CRISPRa kit - CRISPR gene activation of mouse ADP-ribosylation factor-like 4D
Symbol:	Arl4d
Components:	GA211169G1, Arl4d gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211169G2, Arl4d gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211169G3, Arl4d gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q99PE9</u>
Synonyms:	1110036H21Rik; Arf4l; Arfl4; Arl5; ARL6; AW456149
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)