

Product datasheet for **GA211203**

Mouse Tmem108 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem108 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 108
Symbol:	Tmem108
Components:	GA211203G1, Tmem108 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211203G2, Tmem108 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211203G3, Tmem108 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BHE4</u>
Synonyms:	AI462967; B130017P16Rik; R74726
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)