

Product datasheet for GA211206

Mouse Bambi-ps1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Bambi-ps1 CRISPRa kit – CRISPR gene activation of mouse BMP and activin membrane-bound inhibitor, pseudogene (<i>Xenopus laevis</i>)
Symbol:	Bambi-ps1
Components:	GA211206G1, Bambi-ps1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211206G2, Bambi-ps1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211206G3, Bambi-ps1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

