

Product datasheet for GA211218

Mouse Il23a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Il23a CRISPRa kit - CRISPR gene activation of mouse interleukin 23, alpha subunit p19
Symbol:	Il23a
UniProt ID:	<u>Q9EQ14</u>
Synonyms:	IL-23; p19
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)