

Product datasheet for **GA211222**

Mouse **Plekha3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Plekha3 CRISPRa kit - CRISPR gene activation of mouse pleckstrin homology domain-containing, family A (phosphoinositide binding specific) member 3
Symbol:	Plekha3
UniProt ID:	Q9ERS4
Synonyms:	FAPP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)