


Product datasheet for **GA211223**

Mouse Plekha2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plekha2 CRISPRa kit - CRISPR gene activation of mouse pleckstrin homology domain-containing, family A (phosphoinositide binding specific) member 2
Symbol:	Plekha2
Components:	GA211223G1, Plekha2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211223G2, Plekha2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211223G3, Plekha2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9ERS5</u>
Synonyms:	6430512N22Rik; AI851530; AU041791; TAPP2; W91773
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)